Thanks for purchasing our Solight 1T05 Breath Alcohol Computer, a hand held Digital Alcohol Analyzer. Solight 1T05 gives you a graphic indication of your breath alcohol content and thus your possible concentration of intoxication.

Specification

1. Range:	0.1 - 1.5‰ BAC
2. Power Input:	3 x 1.5 AAA Battery
3. Display:	LCD display with bi-color backlighting
4. Auto Power Off:	2 minutes
5. Recycle Time:	~30 seconds
6. Dimension:	107 x 45 x 25 mm
7. Low alcohol level calibration	0.1‰ BAC
8. Accessory:	3 pcs of blowing tubes

1. Install 3 x 1.5V alkaline batteries (AAA) into the battery compartment (auto-cleaning will automatically activate according to the CAUTION manual) The "AUTO CLEANING" indicator flashes and the display counts down. Pressing ON/OFF Key can terminate this "AUTO CLEANING" process.

2. Wash your mouth with water before using. Abstain from smoking or drinking alcohol for at least 5 minutes. Warning!! Using the digital alcohol analyzer immediately or within 5 minutes after consuming alcohol can damage the sensor and distort future measuring. The device is excluded from the guarantee for improper handling.

3. Pressing the ON/OFF Key to turn on the power, sensor pre-heating will be activated for an approximate 9 seconds. "Wait" flashes for 6 seconds and then "PH" & "ST" shows alternatively and the display counts down.

4. Insert a drinking straw or blow tube into the gas inlet located on top of the LCD screen.

5. "READY" lights with a beep tone when count down to zero, and ready to start the test.

6. Press & hold the BLOW Key, the "BLOW" indicator shows on the display. User should blow continuously through the straw for approximately 7 seconds. Stop blowing when you hear a single beep tone and a reading will be displayed on the display.


Alert Signals

1. Low level – readings of less than 0.1‰; the display shows "0.00".
2. Medium level – readings of more than 0.1‰; the analyser shows the alcohol concentration on the display.
3. High level – higher readings; the analyser emits an audible signal, the display goes red, the "NO DRIVING" indicator flashes together with the reading for approx.. 10 seconds.
4. If the reading is more than 1.5‰, the "NO DRIVING" indicator will flash, together with "Hi" (High concentration) with the display lit in red.
5. Press START to turn off the audible signal and start a fresh measurement.

YOU ARE PROHIBITED TO DRIVE since you already exceed the high level.

Notes:

1. To ensure high precision readings, user should perform the auto-cleaning if unused for more than 2 weeks (especially for the first time of using). User should start the auto-cleaning process by removing the batteries and

power up again. The "AUTO CLEANING" indicator flashes. The 5 minutes auto cleaning starts and the display counts down. The 1T05 will give few beep sounds at the last minute of auto cleaning. The analyzer will then turn off automatically after 1 minute and the cleaning process finishes. Before formal test, make sure no alcohol residues left by blowing fresh air into the gas inlet.

2. Please notice that pressing the ON/OFF key at any time can easily terminate the auto-cleaning process.
3. To achieve an accurate reading, user should ensure there is no alcohol residue inside the chamber before starting a new alcohol test. Perform a blank measurement. If the display shows "0.00", everything is working properly. Otherwise, briefly press the START button and then press and hold the BLOW button and blow clean air into the device, or move it backwards and forwards quickly.
4. After any 4 consecutive tests, allow the device to rest for 2 minutes.
5. An "ERROR" flashes on the display for 5 seconds for the following improper operations:
 - When user does not hold the BLOW key for at least 5 seconds during blowing until a beep sound is heard.
 - When user fails to start the test within 5 seconds after "READY" flashes.
6. When the product is unused, please store in a dry environment.
7. When low battery indicator flashes, the reading is not valid, please replace the batteries immediately.
8. If the difference between two successive reading $>0.1\%$ BAC, user should test again.
9. 0.5% BAC = 0.05% BAC = 0.25mg/l (if the convection factor is 2.0)

Many factors (e.g. environment gas residue, smoke & various smell in human mouth, etc...) may affect the BAC level, so the readings are for reference only.

Neither the reseller nor manufacturer accept any liability arising from the use of the unit.

Warning!! Immediate use of the alcohol analyzer or 5 minutes after alcohol consumption can lead to damaging the detection sensor, danger of inaccurate further measurements. The guarantee shall no longer be effective for the device due to inexpert handling

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