### Technické parametry

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### Technical parameters

- measuring range: 32 - 42.9 °C for human body; 0 – 100 °C for object
- distance from the measured object: 3 - 5cm (no fixed distance is necessary)
- working temperature: 10 - 40 ° C - memory of last measurements: 32
- measurement accuracy:  $\leq \pm 0.2$  ° C
- units: °C / °F
- measurement speed: <1 second
- up to 10,000 measurements per pair of batteries (according to quality)
- power supply: 2 pcs AAA (not included in the package)
- automatic shutdown: 15 s

### Battery installation

Observe the correct polarity of the batteries during installation; if you insert the batteries incorrectly, the device will not work.

### Correct use

Correct use of the method is the key to the accuracy of the test, otherwise it may cause measurement errors. Because infrared measurement has high requirements for the surrounding environment, please use it according to the following operations.

1. When measuring the body temperature, the instrument should be pointed to the center of the forehead, above the center of the eyebrow and kept vertical. The measuring part should not be covered by hair. The distance between the instrument and the forehead is recommended to be about 3cm-5cm
2. When the measured person comes from the place with large difference in the measured ambient temperature, he shall stay in the measurement environment for at least 5 minutes, and then measure after the temperature is consistent with the ambient temperature, otherwise the measurement results will be affected.
3. If the forehead of feverish patients is cold, sweating, and other cooling measures are taken, the measurement results will be low, so it should be avoided in this case.
4. When the product is taken out from the place with great difference from the environment to be tested for use, the instrument shall be placed under use for 20 minutes before use.
5. The surrounding environment of the tested person shall be stable, and it is not allowed to test in places with large air flow, such as fan, air conditioner outlet, etc.
6. Do not use the instrument outdoors or in a sunny place.

### Setting

A short press on the "set" button toggles between the "surface" inanimate body measurement mode and the "human body" body temperature measurement mode. The product is specially designed for measuring body temperature, for this reason this mode is preferred. However, the device can also measure, for example, the temperature of food, solids or liquids. Use the "surface" inanimate surface measurement mode for this measurement.

1. Setting the temperature units. Press and hold the "set" button for 5 seconds, the temperature units (Celsius / Fahrenheit) will change on the screen.
2. Review of recent measurements. After each measurement, the device automatically records data, up to the last 32 measurements can be viewed. To view the stored values, select +, -.

### **Presentation of measurement results - measurement of human body temperature**

When the temperature is lower than 37.5 °C, the green backlight of the display will appear (Fig. 1). At a temperature between 37.5 °C - 39 °C, the yellow backlight of the display will light up (Fig. 3). At temperatures above 39 °C, the red backlight of the display appears (Fig. 2).

1

2

3

### **Battery replacement (not included)**

The device can be used up to approx. 10,000 times per battery (depending on the quality). If the battery symbol flashes on the display, it means that the battery is low and needs to be replaced. Do not use rechargeable batteries or use this product with one battery for more than half a year

### **Routine maintenance**

1. Pay attention to external cleaning, and wipe the dirt with alcohol or disinfectant;
2. Internal dirt: do not wipe with hands, wipe the lens surface with cotton swab stained with more than 95% anhydrous alcohol;

Storage: please keep in a dry and dark place out of direct sunlight

### **Troubleshooting**

1. The temperature displayed on the screen is greater than 95 °C, and the temperature is Fahrenheit, which is converted into Celsius
2. The temperature mode of the screen display itself is lower than 32 degrees, and the temperature mode should be converted into a body in the "surface temperature" mode
3. „Hi“ information is displayed on the screen. In the temperature mode, when the measured temperature is higher than 42.9 degrees, the product will display hi. First, confirm the working mode of the product. In the temperature mode, the measured human body will still display hi. First, check whether the measured part of the human body has contacted the external heat source. If not, contact our technical personnel
4. The screen displays „Lo“ information. In body temperature mode, when the measured temperature is lower than 32 °C, the instrument will display lo. Confirm that the instrument is in surface temperature or body temperature mode, and the measured human body still displays lo. First, find out the reason according to the following table. If it is confirmed that it is not the reason in the table, contact the technical personnel.

The main reasons for the display of „Lo“ are as follows:

- There is hair or sweat when reading the temperature;
- There is cold air blowing on the forehead;
- Make sure there is no obstacle or sweat on the forehead;
- The forehead has just been cold applied;
- The measuring distance is too far

The product is issued with a CE declaration of conformity in accordance with applicable regulations. On request from the manufacturer: [info@solight.cz](mailto:info@solight.cz), or for download at [www.solight.cz](http://www.solight.cz). Manufacturer: Solight Holding, s.r.o., Na Brně 1972, Hradec Králové 500 06, Czech Republic

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